



Natasha wears
two Rose Gold
MCPDrivers.
[right]

Midnight Moss
PIPDriver.
[below]



npdevices.com

NEW METALLIC COATINGS FOR ALL DEVICES

Introducing **Rose Gold** and **Midnight Moss** coatings for the PIPDriver, MCPDriver, ThumbDriver, and GripLock Finger.



PIPDriver™

Custom device driven by an intact PIP joint with enough residuum to engage the ring.



MCPDriver™

Custom device driven by an intact MCP joint with enough residuum to engage the ring.



ThumbDriver™

Custom device driven by intact MCP and CMC joints in the thumb. However, patients with more proximal amputations may still benefit from the device as long as their residuum can engage the ring.

GripLock Finger™

Elegant solution for amputations proximal to the MCP joint. Positioned using a surface or their opposing hand, and then return to full extension by releasing the dorsal pawl or moving to the full flexion release point.

