

## PIPDriver | 2nd Generation

### OTHER COMFORT ASPECTS INCLUDE:

- Custom designed to fit comfortably and suit your patient's essential daily needs.
- Easy to clean, care for, don and doff.
- Made from strong and durable medical-grade nylon 12.
- Cage-like structure protects residuum.
- Slender and sleek design allows for wearing of multiple devices.
- Conforming silicone rubber tip pads resemble the natural fingertip.
- Lightweight titanium fasteners.



### Function is what the PIPDriver is all about.

Designed for individuals amputated at the middle phalanx. Of the products we offer, the **PIPDriver** is the most simple and intuitive to use. Because the device fits alongside the patient's finger, we can anatomically match their PIP and DIP joints so it operates predictably and naturally. The **PIPDriver** offers exceptional daily utility.

The second generation **PIPDriver** combines the same reliable functionality with a new, sleek aesthetic design, improved hinges, and a variety of color options.

### IDEAL CANDIDATE SPECIFICATIONS:


- This device is driven by an intact PIP joint with enough residuum to engage the ring.

For an amputee who doesn't meet the minimum requirements for the **PIPDriver**, it is possible to obtain some functional benefit from the **PIPDriver**. However, special considerations will have to be made prior to ordering a device.

Sizing is as simple as ordering a Naked Prosthetics sizing kit, completing the forms, submitting detailed pictures of both hands and a video showing range of motion. From this information, our engineers design and manufacture a custom-functional, finger prosthesis.

The **PIPDriver** has a four to six week lead time from the point all required documentation is received.

### AVAILABLE COATINGS:



Organic



Minimalist